

Run-to-Run (Indiv. & MR Charts)

General Results

- Stable within each Trial.
- Trial 1 to Trial 2, Merit & Next - product shift, new batch.
- **Out of Control** -
 - All machines, Trial 1 vs. Trial 2, Next, Nicotine. Product batch difference.
 - M59, Trial 1 vs. Trial 2, Next, TAR. Product batch difference?
 - M59, Run 17, Merit, high Nicotine.
See S chart results (variation high), probably product or GC.
 - M58, Run 16, Monitor, low H₂O & high Tar. Lab? Holiday?
 - M59, Runs 4 & 5, Monitor, low TPM & Tar, high Nicotine. Product variation?
- Batch differences not as pronounced in Tar.